

**PRECAUTIONARY STATEMENTS:
HAZARDS TO HUMANS AND
DOMESTIC ANIMALS**

Caution

Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

FIRST AID

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after first 5 minutes, then continue rinsing eye.

Call a poison control center or doctor for treatment advice.

Have the product container or label with you when calling poison control center or doctor or going for treatment.

Use trophon® Companion Cleaning Wipes in:

- Hospitals
- Medical offices
- Nursing homes
- Health clubs
- Schools
- Ambulances
- Emergency vehicles

Use trophon Companion Cleaning Wipes on washable hard, nonporous surfaces of:

- Computer keyboards
- Countertops
- Desks
- Tiles, glazed
- Walls

Other hard nonporous surfaces made of: Glass, Glazed ceramic, Glazed porcelain, LDPE, HDPE, Metal, Plastic, PVC, Stainless steel, Upholstery – vinyl and plastic.

Non-critical medical device surfaces:

- Bed frames
- Bed rails
- Beds
- Bladder scan
- Dopplers
- Medical examining tables
- Non-critical hospital equipment surfaces
- Radiology equipment
- Stethoscopes
- Ultrasound monitors
- Ultrasound transducers and equipment

trophon Companion Cleaning Wipes kills:

Bacteria:

- Pseudomonas aeruginosa
- Salmonella enterica
- Staphylococcus aureus
- Acinetobacter baumannii
- Burkholderia cepacia
- Campylobacter jejuni
- Corynebacterium ammoniagenes
- Enterobacter aerogenes
- Enterococcus faecalis
- Enterococcus faecalis - Vancomycin Resistant (VRE)
- Escherichia coli
- Escherichia coli O157:H7
- ESBL Escherichia coli
- Klebsiella pneumoniae
- Klebsiella pneumoniae (NDM-1)
- Legionella pneumophila
- Listeria monocytogenes
- Salmonella schottmuelleri
- Salmonella typhi
- Serratia marcescens
- Shigella dysenteriae
- Staphylococcus aureus - Methicillin-Resistant (MRSA)
- Staphylococcus aureus - Multi-Drug Resistant
- Staphylococcus aureus - Vancomycin Intermediate Resistant - (VISA)
- Streptococcus pyogenes
- Vibrio cholerae

Viruses:

- *Hepatitis B Virus
- *Hepatitis C Virus
- *Herpes Simplex Virus Type 1
- *Herpes Simplex Virus Type 2
- *HIV-1
- *Influenza Type A Virus / Brazil
- *Influenza A Virus H1N1
- *2013 Influenza A virus (H7N9)
- *Norwalk Virus (Norovirus)
- *Respiratory Syncytial Virus
- *Rotavirus
- *SARS Associated Coronavirus (SARS)
- *SARS-Related Coronavirus 2, cause of COVID-19
- *Vaccinia Virus

Animal Viruses:

- Avian Influenza Virus H3N2
- Avian Influenza Virus H5N1
- Canine Distemper Virus
- Feline Calicivirus
- Newcastle's Disease Virus
- Pseudorabies Virus

NET CONTENTS:
• 100 wipes
• 11 in x 7 in
• 1.92 lbs

trophon®

**Companion
Cleaning Wipes**

Cleans in one easy step

- No rinsing required
- Multi-purpose
- Multi-surface
- Convenient one-hand dispensing
- Everyday Disinfectant Cleaner
- Fast drying; leaves no dull or sticky residue
- Kills 99.9% of Bacteria and Viruses*
- This product will clean and disinfect hard, non-porous, non-food surfaces.

Active Ingredients:

Octyl decyl dimethyl ammonium chloride.....	0.0909%
Dioctyl dimethyl ammonium chloride.....	0.0364%
Didecyl dimethyl ammonium chloride.....	0.0545%
Alkyl (C14, 50%; C12, 40%; C16, 10%) dimethyl benzyl ammonium chloride.....	0.1212%
Other Ingredients:	99.6970%
Total.....	100.0000% ¹

¹Includes weight of absorbent wipe.

**KEEP OUT OF
REACH OF
CHILDREN
CAUTION**



DO NOT FLUSH

See left panel for Precautionary Statements and First Aid

nanosonics

To Open Package:

Remove cover. Pull up wipe corner from center of roll. Thread through hole in the cover. Replace cover. Open center cap on cover to remove wipe. When not in use, close center cap securely to prevent moisture loss.

DIRECTIONS FOR USE:

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Do not use as a diaper wipe or for personal cleansing. This is not a baby wipe! Do not flush in toilet.

For plastic and painted surfaces, spot test on an inconspicuous area before use. **trophon Companion Cleaning Wipes** will not harm most surfaces: including acrylic, sealed fiberglass and vinyl.

TO CLEAN AND DEODORIZE:

Wipe surface with wipe. Let air dry. Discard used wipe in trash. Do not flush in toilet.

CLEANING AND DISINFECTION/VIRUCIDAL* DIRECTIONS:

Wipe hard, nonporous surface with wipe until surface is visibly wet. Use enough wipes to keep surfaces visibly wet for 4 minutes.

Allow surface to remain visibly wet for 4 minutes.

For Norwalk Virus, Feline Calicivirus, Norovirus and Rotavirus, treated surfaces must remain visibly wet for 10 minutes. Let air dry.

If surfaces are visibly dirty, clean first with another wipe before disinfecting. Discard used wipe in trash. Do not flush in toilet.

***KILLS HIV AND HBV AND HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES / OBJECTS PREVIOUSLY SOILED WITH BLOOD / BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces / objects with blood or body fluids, and in which the surfaces / objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 [HIV-1] [associated with AIDS] or Hepatitis B Virus [HBV] or Hepatitis C Virus [HCV].

Special Instructions for Cleaning and Decontamination against HIV-1 or HBV or HCV on Surfaces/Objects Soiled with Blood/Body Fluids

Personal protection:

Clean-up must always be done wearing protective gloves, gowns, masks and eye protection.

Cleaning procedure:

Blood and other body fluids containing HIV or HBV or HCV must be thoroughly cleaned from surfaces and objects before application of **trophon® Companion Cleaning Wipes**. **trophon® Companion Cleaning Wipes** can be used for this purpose.

Contact time:

Leave surface visibly wet for 1 minute for HIV-1. Use appropriate contact time as listed on label for all other viruses and bacteria claimed.

Disposal of infectious material:

Blood, body fluids, cleaning materials and clothing must be disposed of according to local regulations for infectious waste disposal.

This product is for use on hard nonporous surfaces not conducive to treatment by immersion or excess liquid against Avian Influenza A.

This product is not to be used as a terminal sterilant / high level disinfectant on any surface or instrument that is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

STORAGE AND DISPOSAL

Store in original container in areas inaccessible to children. Keep securely closed. Nonrefillable container. Do not reuse or refill this container. Discard empty container in trash or recycle if available.

EPA Reg. No. 6836-340-93267

EPA Est. No. A=74058-CA-001 B=74058-OH-001

Nanosonics, Inc.
7205 E. 87th St, Indianapolis, Indiana 46256
1-844-876-7466
customerservice@nanosonics.us
www.nanosonics.us

Manufactured for : Nanosonics Limited
7-11 Talavera Road
Macquarie Park
NSW 2113, Australia
+61 2 8063 1600
customerservice@nanosonics.com.au
www.nanosonics.com.au

REF N00112-US

L00140 6.0