

#### To open package:

Remove cover. Pull up wipe corner from center of roll. Thread through hole in the cover. Replace cover. Open center cap on cover to remove wipe. When not in use, close center cap securely to prevent moisture loss.

#### DIRECTIONS FOR USE:

It is a violation of federal law to use this product in a manner inconsistent with its labeling.

Do not use as a diaper wipe or for personal cleansing.

This is not a baby wipe!

Do not flush in toilet.

For plastic and painted surfaces, spot test on an inconspicuous area before use.

Do not use on unpainted wood.

#### TO CLEAN:

Wipe surface with wipe. Let air dry. Discard used wipe in trash. Do not flush in toilet.

#### CLEANING AND DISINFECTION/VIRUCIDAL\*

##### DIRECTIONS:

Wipe surface with wipe until surface is visibly wet.

Use enough wipes to keep surfaces wet for 4 minutes.

For *Candida albicans*, treated surfaces must remain wet for 8 minutes.

For norwalk virus, feline calicivirus, norovirus and rotavirus, treated surfaces must remain wet for 10 minutes.

Let air dry. If surfaces are extremely dirty, clean first with another wipe before disinfecting. Discard used wipe in trash. Do not flush in toilet.

Rinse food contact surfaces with clean, potable water after disinfecting. Do not use to disinfect dishes, glassware or utensils.

#### \*KILLS HIV AND HBV AND HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES / OBJECTS PREVIOUSLY SOILED WITH BLOOD / BODY FLUIDS

in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces / objects with blood or body fluids, and in which the surfaces / objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus type 1 [HIV-1] [associated with AIDS] or hepatitis B virus [HBV] or hepatitis C virus [HCV].

**Special instructions for cleaning and decontamination against HIV-1 or HBV or HCV on surfaces / objects soiled with blood / body fluids.**

**Personal protection:** Clean-up must always be done wearing protective gloves, gowns, masks and eye protection.

**Cleaning procedure:** Blood and other body fluids containing HIV or HBV or HCV must be thoroughly cleaned from surfaces and objects before application of **trophon® Companion Cleaning Wipes**. **trophon® Companion Cleaning Wipes** can be used for this purpose.

**Contact time:** Leave surface wet for 1 minute for HIV-1.

Use appropriate contact time as listed on label for all other viruses, bacteria and fungi claimed.

**Disposal of infectious material:** Blood, body fluids, cleaning materials and clothing must be disposed of according to local regulations for infectious waste disposal.

This product is for use on hard nonporous surfaces not conducive to treatment by immersion or excess liquid against avian influenza A.

This product is not to be used as a terminal sterilant / high level disinfectant on any surface or instrument that is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

Use **trophon® Companion Cleaning Wipes** in:

Hospitals Schools  
Medical offices Nursing homes  
Ambulances Emergency vehicles  
Health clubs

Use **trophon® Companion Cleaning Wipes** on washable hard, nonporous surfaces of:  
Computer keyboards • Tiles, glazed • Countertops  
• Walls • Desks

*Non-critical medical device surfaces:*

Beds  
Bed frames  
Bladder scan  
Bed rails  
Dopplers  
Non-critical hospital equipment surfaces  
Medical examining tables  
Radiology equipment  
Stethoscopes  
Ultrasound monitors  
Ultrasound transducers and equipment

*Other hard nonporous surfaces made of:*

Glass  
Glazed ceramic  
Glazed porcelain  
LDPE, HDPE  
Metal  
Plastic  
PVC  
Stainless steel  
Upholstery- vinyl and plastic  
Note: Do not use on unpainted wood or natural marble.

# trophon®

## Companion Cleaning Wipes

### Cleaner – Disinfectant

#### One step cleaning • Larger-sized wipe • No rinsing required

- Fragrance-free
- Fast drying; leaves no dull or sticky residue
- Multi-purpose
- Multi-surface
- Non-abrasive
- This product will clean and disinfect hard, non-porous surfaces.
- Kills 99.9% of bacteria and viruses\*
- Convenient one-hand dispensing
- Recyclable packaging

Active Ingredients:

Octyl decyl dimethyl ammonium chloride.....	0.0909%
Diocetyl dimethyl ammonium chloride.....	0.0364%
Didecyl dimethyl ammonium chloride.....	0.0545%
Alkyl (C <sub>14</sub> , 50%; C <sub>12</sub> , 40%; C <sub>16</sub> , 10%) dimethyl benzyl ammonium chloride.....	0.1212%
Other Ingredients:.....	99.6970%
Total.....	100.0000% <sup>1</sup>

<sup>1</sup> Includes weight of absorbent wipe.

## KEEP OUT OF REACH OF CHILDREN

# CAUTION

See right panel for precautionary statements and first aid

Net Contents: 100 wipes  
11" x 7" (27.94cm x 17.78cm)  
2lb 7oz (1106g)

EPA Reg. No. 6836-340-93267  
EPA Est. No. 73027-WI-1



#### trophon® Companion Cleaning Wipes kills:

Bacteria:

*Pseudomonas aeruginosa*  
*Salmonella enterica*  
*Staphylococcus aureus*  
*Acinetobacter baumannii*  
*Burkholderia cepacia*  
*Campylobacter jejuni*  
*Corynebacterium ammoniagenes*  
*Enterobacter aerogenes*  
*Enterococcus faecalis*  
*Enterococcus faecalis*- vancomycin resistant (VRE)  
*Escherichia coli*  
*Escherichia coli* O157:H7  
ESBL *Escherichia coli*  
*Klebsiella pneumoniae*  
*Klebsiella pneumoniae* (NDM-1)  
*Legionella pneumophila*  
*Listeria monocytogenes*  
*Salmonella schottmuelleri*  
*Salmonella typhi*  
*Serratia marcescens*  
*Shigella dysenteriae*  
*Staphylococcus aureus* - methicillin-resistant (MRSA)  
*Staphylococcus aureus* - multi-drug resistant  
*Staphylococcus aureus* - vancomycin intermediate resistant – (VISA)  
*Streptococcus pyogenes*  
*Vibrio cholerae*

Viruses:

\*Hepatitis B virus  
\*Hepatitis C virus  
\*Herpes simplex virus type 1  
\*Herpes simplex virus type 2  
\*HIV-1  
\*Human coronavirus  
\*Influenza type A virus / Brazil  
\*Influenza A virus H1N1  
\*2013 influenza A virus (H7N9)  
\*Norwalk virus (norovirus)  
\*Respiratory syncytial virus  
\*Rotavirus  
\*SARS associated coronavirus (SARS)  
\*Vaccinia virus (pox virus)

Animal viruses:

Avian influenza virus H3N2  
Avian influenza virus H5N1  
Canine distemper virus  
Feline calicivirus  
Newcastle's disease virus  
Pseudorabies virus

Fungi:

*Candida albicans*

#### PRECAUTIONARY STATEMENTS

#### HAZARDS TO HUMANS AND DOMESTIC ANIMALS

Caution. Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

#### FIRST AID

**IF IN EYES:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after first 5 minutes, then continue rinsing eye.

Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling poison control center or doctor or going for treatment.

## STORAGE AND DISPOSAL

Store in original container in areas inaccessible to children. Keep securely closed. Nonrefillable container. Do not reuse or refill this container. Discard empty container in trash or recycle if available.

Manufactured for:  
Nanosonics Limited  
14 Mars Road, Lane Cove  
NSW 2066, Australia  
+61 2 8063 1600  
info@nanosonics.com.au

USA contact:  
Nanosonics, Inc.  
7205 E. 87th Street  
Indianapolis, IN 46256, USA  
1-844-876-7466  
support@nanosonics.us  
www.nanosonics.us

REF N00112-US

