trophon®

Wireless Ultrasound Probe Holder (WUPH-A)

Nanosonics Limited
7-11 Talavera Road
Macquarie Park
NSW 2113, Australia
+61 2 8063 1600
customerservice@nanosonics.com.au
www.nanosonics.com.au

Nanosonics, Inc.
7205 E. 87th Street
Indianapolis, Indiana 46256
1-844-876-7466
customerservice@nanosonics.us
www.nanosonics.us


Made in China

REF N05030-US

nanosonics
Infection Prevention. For Life.

L03258 3.0 10/2024

Figure 1

trophon® Wireless Ultrasound Probe Holder

trophon® AcuTrace® Medical Instrument Tag

Figure 2a


Figure 2b

Figure 2c

Figure 3


Figure 4a - trophon2


Figure 4b - trophon EPR


Figure 5


Figure 6



Figure 7

ENGLISH

⚠ Only the wireless probes shown in Table 1 are compatible with trophon® Wireless Ultrasound Probe Holder.

trophon Wireless Ultrasound Probe Holder is hereafter referred to as the WUPH-A in the rest of this document.

Table 1: WUPH-A Compatible Probes List

Manufacturer	Probe Type	Model
GE Healthcare	Dual Window	Vscan Air™ CL[1] / SL[1]
Siemens	Linear	ACUSON Freestyle L8-3[2] / L13-5[2] / L17-5[2]
Siemens	Curved	ACUSON Freestyle C5-2[2,3] / CH5-2[2,3]

Table 1 Legend:
[1] Probe with dual windows: Ensure the wider end is facing upwards after inserting it into the WUPH-A. See Figure 2a.

[2] Probe with the removable battery: Ensure the battery is removed and load the probe with the window end facing down. See Figure 2b.

[3] Curved probe with wider window: Load the probe with the window facing down at an angle. After the window passes through, straighten the probe so the probe body is resting on all of the contact points. See Figure 2c.

Intended Use: trophon® Wireless Ultrasound Probe Holder, when used in accordance with its labelling, is a non-sterile wireless ultrasound probe holder. It is intended to hold a wireless ultrasound probe in the trophon chamber.

Indications for Use:

The trophon Wireless Ultrasound Probe Holder is an accessory to be added to the trophon EPR and trophon2 to allow wireless ultrasound probes to be disinfected.

trophon EPR:
The trophon EPR is designed to provide High-Level Disinfection of ultrasound transducers. The system uses the trophon Disinfectant which is intended to be used exclusively with the trophon EPR device.

The trophon Disinfectant is intended for use as a High-Level Disinfectant to be used exclusively with the trophon EPR for the High-Level Disinfection of ultrasound transducers.

The trophon EPR is suitable for use in general hospital and health care facilities by trained personnel.

The trophon EPR system consists of a multiple use instrument combined with a single use disinfectant, delivered from a multi-dose cartridge.

The trophon Disinfectant should be used with the following contact conditions:
Minimum Operational Cycle Time: 7 minutes
Minimum Concentration: 31.5%
Minimum Disinfectant Dose: 1.0 g
Minimum Chamber Temperature: 56°C

trophon2:
The trophon2 is designed to provide High-Level Disinfection (HLD) of validated ultrasound transducers. High-Level Disinfection is achieved by surface exposure to a controlled dose of hydrogen peroxide mist delivered to a disinfection chamber contain the ultrasound probe.

The trophon2 system consists of a multiple use instrument combined with a single use disinfectant "trophon Sonex-HL", delivered from a multi-dose cartridge.

The trophon2 is suitable for use in general hospital and health care facilities by trained personnel.

The trophon Sonex-HL should be used with the following contact conditions:
Minimum Operational Cycle Time: 7 minutes
Minimum Concentration: 31.5%
Minimum Disinfectant Dose: 1.0 g
Minimum Chamber Temperature: 56°C

Figure 1 Legend:
1. Top loop
2. Handle
3. Contact points

Visual Inspection of the WUPH-A: The WUPH-A shall be visually inspected prior to use to ensure it is clean, dry, and not damaged. Damage to the WUPH-A includes changes to surface smoothness or shape, and appearance of cracks, holes and chips shown in Figure 5. If the WUPH-A is damaged, follow the 'Disposal' instructions below.

Cleaning the WUPH-A: Wipes containing quaternary ammonium compounds have been validated for cleaning efficacy and chemical compatibility of the WUPH-A, and can be used to clean the WUPH-A. Clean the WUPH-A and ensure the wipe contacts all surfaces (outer, inner surfaces and handle). Other cleaning products may be used. Please consult the IFU of cleaning products for suitability. Ensure that the WUPH-A has been thoroughly dried after cleaning with a dry, clean, single-use, low-lint, or lint-free wipe, such as the trophon Companion Drying Wipes (**REF** N00113) and ensure that no visible debris is present. The cleaning process should be repeated if the WUPH-A is still not visually clean. If the WUPH-A cannot be cleaned, it should be safely disposed following the 'Disposal' instructions below.

Disinfecting the WUPH-A: trophon devices have been validated for high-level disinfection of the WUPH-A. The WUPH-A is disinfected during use in the trophon device.

Recommended Product Lifetime: The WUPH-A has been verified to last a minimum of 5,000 trophon cycles when adequately handled and reprocessed. The WUPH-A should be visually inspected for damage prior to use following the 'Visual Inspection of the WUPH-A' instructions.

trophon AcuTrace Medical Instrument Tag (trophon2 only)
Refer to the trophon2 User Manual for detailed instructions on the trophon AcuTrace Medical Instrument Tag digital traceability system.

The trophon AcuTrace Medical Instrument Tag can be programmed to store a combination of the probe and WUPH-A identification information using the trophon2 touchscreen. To implement this, it is recommended that the input for the Probe Name includes both the probe and the WUPH-A information such as the last digits of each Unique Device Identifier (UDI).

The trophon AcuTrace Medical Instrument Tag can be attached to the top loop of the WUPH-A. Up to three tags can be added to one WUPH-A. Refer to the trophon2 User Manual for detailed instructions on digital traceability with AcuTrace technology.

Setup and Use
Refer to the trophon EPR or trophon2 User Manuals for detailed instructions on how to use the trophon device.

⚠ Only the wireless probes shown in Table 1 are compatible with WUPH-A.

1. Wear clean gloves.
2. Visually inspect the WUPH-A to ensure it is clean, dry, and not damaged following the 'Visual Inspection of the WUPH-A' instructions.

⚠ Do not use the WUPH-A if it is visually damaged (as shown in Figure 5). If the WUPH-A appears to be visually damaged, follow the 'Disposal' instructions.

⚠ Do not use the WUPH-A if it is dirty or wet (as shown in Figure 6). If cleaning is required, follow the 'Cleaning the WUPH-A' instructions.

⚠ Do not attach anything to any part of the WUPH-A other than the handle top loop, as shown in Figure 7.

3. Clean and dry the probe and check carefully for probe defects prior to commencing the high-level disinfection (HLD) process, as per the probe manufacturer's instructions. Ensure that the probe has been thoroughly dried after cleaning with a dry, clean, single-use, low-lint or lint-free wipe and that no visible debris is present.
4. Place the probe in the WUPH-A: With one hand, hold the probe with the wider end facing upwards and using the other hand, insert the probe into the holder, from the top as shown in Figure 2a, 2b and 2c.

Ensure the probe is resting on all four touch points of the WUPH-A as shown in Figure 3.

5. If AcuTrace is enabled, when prompted on the trophon2 screen, scan the trophon AcuTrace Medical Instrument Tag that corresponds to the associated probe against the AcuTrace reader.
6. Open the trophon chamber door and insert the WUPH-A containing the probe into the trophon EPR/trophon2 device, as shown in Figure 4. Ensure the handle of the WUPH-A is appropriately positioned in the gland seal with the trophon logo facing outwards. Refer to Figures 4a and 4b for the final position in trophon2 and trophon EPR devices, respectively.

NOTE: The WUPH-A may sit lower than shown in Figure 4b depending on the type of trophon EPR clamp used.

⚠ Incorrect positioning of the probe in the WUPH-A or incorrect positioning of the WUPH-A in the trophon device may result in unsuccessful HLD cycles.

7. Follow the trophon EPR or trophon2 User Manual for 'Inserting the Chemical Indicator', 'Closing the Chamber Door' and 'Disinfecting the Probe'.
8. Replace gloves.
9. When the cycle is complete, carefully remove the WUPH-A from the trophon device with two hands. Handle with care to avoid damage and pay attention not to drop the probe.
10. Wipe the probe with a dry, clean, single-use, low-lint or lint-free wipe. Visually inspect the probe and remove any present disinfectant residue as shown in Figure 6.
11. Remove the used trophon Chemical Indicator and verify the color change against the Color Assessment Chart on the trophon Chemical Indicator carton. Discard the used trophon Chemical Indicator.
12. Refer to the trophon2 User Manual for 'Confirming the HLD cycle'.
13. If the trophon EPR/trophon2 device indicates a passed cycle, the probe is now ready for use or storage in a Clean Ultrasound Probe Cover (**REF** N00102).
14. The WUPH-A should be stored in a clean and dry area away from potential contamination or damage. To reduce the risk of contamination, the WUPH-A can be stored in a Clean Ultrasound Probe Cover (**REF** N00102).

Do not store the WUPH-A when it is wet. When required, wipe the WUPH-A with a dry, clean, single-use, low-lint or lint-free wipe.

Disposal: Consult with your local environmental and government agencies to confirm how the WUPH-A should be disposed, such as into general waste or recycling, and any local rules that apply in your area.

NOTE: The WUPH-A is not compatible with the trophon Curved Probe Positioner (CPP).

⚠ Any serious incidents in relation to trophon Wireless Ultrasound Probe Holder should be reported to Nanosonics and/or the authority of your member state.

Glossary

	Manufacturer	REF	Product Number
EC REP	European Authorized Representative	SN	Serial Number
	EU Importer		Date of Manufacture
LOT	Batch No.	CE	Conforming to Medical Devices Regulation (EU)2017/745 (MDR)
	Consult Instructions for Use	MD	Medical Device
	Waste Electrical and Electronic Equipment		Restriction of Hazardous Substances
UK CA	Confirming to UK MDR		

trophon®
Wireless Ultrasound Probe Holder (WUPH-A)

Nanosonics Limited
7-11 Talavera Road
Macquarie Park
NSW 2113, Australia
+61 2 8063 1600
customerservice@nanosonics.com.au
www.nanosonics.com.au

Nanosonics, Inc.
7205 E. 87th Street
Indianapolis, Indiana 46256
1-844-876-7466
customerservice@nanosonics.us
www.nanosonics.us


Made in China

REF N05030-US

nanosonics
Infection Prevention. For Life.

L03258 3.0 10/2024

Figure 1

trophon® Wireless Ultrasound Probe Holder

trophon® AcuTrace® Medical Instrument Tag

Figure 2a

Figure 2b

Figure 2c

Figure 3


Figure 4a - trophon2


Figure 4b - trophon EPR


Figure 5


Figure 6



Figure 7

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