

## **trophon Sonex-HL Ingredient Communication**

California SB-258 Cleaning Products Right to Know Act of 2017

Full Disclosure of All Intentionally Added and Nonfunctional Ingredients

<b>trophon Sonex-HL</b>		
<b>Chemical Name</b>	<b>CAS No.</b>	<b>Function</b>
<b>Water</b>	<b>7732-18-5</b>	<b>Diluent</b>
<b>Hydrogen Peroxide</b>	<b>7722-84-1</b>	<b>Active Ingredient</b>
<b>Disodium pyrophosphate</b>	<b>7758-16-9</b>	<b>Stabilizer</b>