

## **trophon Companion Cleaning Wipes Ingredient Communication**

California SB-258 Cleaning Products Right to Know Act of 2017

Full Disclosure of All Intentionally Added and Nonfunctional Ingredients

<b>trophon Companion Cleaning Wipes</b>		
<b>Chemical Name</b>	<b>CAS No.</b>	<b>Function</b>
<b>Water</b>	<b>7732-18-5</b>	<b>Diluent</b>
<b>Benzalkonium Chloride</b>	<b>68424-85-1</b>	<b>Active Ingredient</b>
<b>Quaternium-24</b>	<b>32426-11-2</b>	<b>Active Ingredient</b>
<b>Didecyldimonium Chloride</b>	<b>7173-51-5</b>	<b>Active Ingredient</b>
<b>Tetrasodium EDTA</b>	<b>64-02-8</b>	<b>Chelating Agent</b>
<b>Dimethyldioctylammonium Chloride</b>	<b>5538-94-3</b>	<b>Active Ingredient</b>
<b>Sodium Silicate Pentahydrate</b>	<b>10213-79-3</b>	<b>Detergent Builder</b>
<b>Polyaminopropyl Biguanide</b>	<b>32289-58-0</b>	<b>Preservative</b>