

trophon²

Protection by design



High-level disinfection for ultrasound probes

Consumables and accessories



With AcuTrace[®] (RFID) technology to enhance digital record-keeping and deliver easy, reliable workflows for ultrasound probe high-level disinfection (HLD).

nanosonics
Infection Prevention. For Life.



High-level disinfection device

trophon®2 Device

The latest innovation in ultrasound probe high-level disinfection featuring an enhanced design, streamlined workflows, plus AcuTrace RFID technology for supporting digital record traceability.



Consumables

trophon Sonex HL®

Sonex-HL is a proprietary disinfectant liquid with 35% hydrogen peroxide chemistry to achieve effective high-level disinfection of the entire ultrasound probe including the shaft and handle. Incorporates AcuTrace technology to store batch and expiry details directly into trophon2.



trophon Chemical Indicator

The trophon Chemical Indicator responds to a combination of disinfectant dosage and disinfectant concentration. These parameters need to be met for the chemical indicator to change to a pass color.


This qualitative color change provides additional verification of success of each disinfection cycle. Indicator packaging incorporates AcuTrace technology to record and store batch and expiry details directly onto your trophon2 device.



Traceability accessories

trophon AcuTrace Operator Card

Incorporates AcuTrace technology to digitally capture operator interaction details, stored directly on your trophon2 device.



trophon AcuTrace Medical Instrument Tag

Incorporates AcuTrace technology to digitally capture ultrasound probe details, stored directly on your trophon2 device.



trophon Printer

The trophon Printer is a convenient traceability solution, printing high-quality trophon disinfection labels each and every cycle.



trophon Printer Label Roll

High quality printer labels with a minimum 7-year life span, appropriate for medical records.



trophon Logbook

Simplifies and systematizes record-keeping while providing evidence of quality HLD reprocessing at healthcare facilities.

Point-of-care set-up accessories



trophon Wall Mount 2

A convenient option for installing a trophon2 device where there are space constraints.

trophon Printer Wall Mount[^]

Custom designed for secure, horizontal mounting of the trophon Printer to a wall. Includes a cable tidy area for convenience.



trophon Cart

Makes trophon2 devices fully mobile for convenient point of care use and adaptability where there are space constraints.

trophon Printer Cart Mount 2[^]

Designed for securely attaching the trophon Printer to the trophon Cart, in an easily accessible, convenient location.

trophon Cart UPS Bracket

Designed to accommodate an Uninterruptible Power Supply (UPS) unit to enable trophon devices mobility between patient examination rooms eliminating the need to turn off the trophon device.*

Wipes



trophon Companion Cleaning Wipes

Ideal for use with trophon2 devices for cleaning an ultrasound probe before HLD. The multi-purpose wipes also provide fast and efficient low-level disinfection (LLD) of ultrasound probes and other equipment including consoles, monitors, keyboards, cables and probe holders.



trophon Companion Drying Wipes

Used to dry probes before and after high-level disinfection. trophon Companion Drying Wipes are highly-absorbent to leave surfaces completely dry.

trophon Companion Wipes Wall Bracket

For wall-mounting trophon Companion Drying Wipes for easier access and dispensing.

Ultrasound probe storage cover



Clean Ultrasound Probe Cover

Custom designed probe covers to prevent endocavitary and surface ultrasound probes from recontamination after high level disinfection and before they are used on the next patient.

[^]Printer not included.

*UPS not included. Contact your local sales representative for compatible UPS specifications. Disinfection cycles should not be run while the trophon device is unplugged from the mains power.

trophon2 System Product Descriptions

Product Name	Product Code	Quantity
trophon2	N05000-US	1 Device
trophon Sonex-HL	N05002	1 Box of 6 Cartridges Minimum order quantity: 2 boxes
trophon Chemical Indicator	N05003	1 Box of 300 Chemical Indicators
trophon AcuTrace Operator Card	N05006	1 Box of 10 Cards
trophon AcuTrace Medical Instrument Tag	N05007	1 Box of 10 Tags
trophon Printer	N00048-NNA	1 Printer
trophon Printer Wall Mount	N00105	1 Mount
trophon Printer Label Roll	N00049-ROW	1 Roll of 600 Labels*
trophon Printer USB Cable	N05011	1 Cable
trophon Logbook	N00098	1 Box of 5 Logbooks
trophon Wall Mount 2	N05009	1 Mount
trophon Cart	N05010	1 Box of 1 Cart
trophon Printer Cart Mount 2	N05013	1 Mount
trophon Cart UPS Bracket	N00110	1 Bracket
trophon Companion Cleaning Wipes	N00112-US	1 Box of 6 Canisters x 100 Wipes
trophon Companion Drying Wipes	N00113	1 Box of 6 Canisters x 100 Wipes
trophon Companion Wipes Wall Bracket	N00114	1 Bracket
Clean Ultrasound Probe Covers	N00102	1 Box of 100 Covers

*1 Roll of 350 labels for Parametric Release.



Nanosonics, Inc.
(Distributor USA)
7205 E 87th Street
Indianapolis, IN 46256, USA
T: 1-844-876-7466
E: info@nanosonics.us
W: www.nanosonics.us

Nanosonics Limited
(Manufacturer)
7-11 Talavera Road
Macquarie Park, NSW 2113, Australia
T: +61 2 8063 1600
E: info@nanosonics.com.au
W: www.nanosonics.com.au